

ILUMIASKIN®



Skin care devices market is projected to be valued at USD \$23.8 Billion (AUD \$33.3 Billion) in 2023 and estimated to grow to USD \$101.5 Billion (AUD \$142.1 Billion) in 2033.¹

1. Retail Beauty Autumn 2023. **2.** Allied Market Research. **3.** Fior Markets 2020 Australia LED Face Masks Market -2027.

The global LED mask market size was valued at USD \$267.3 million in 2022, and is projected to reach USD \$702.9 million

by 2032.²






The home-use LED mask/ market is significant and continuing to grow.³ Currently valued at AUD ~ \$9.4 MILLION. Expected to grow to AUD \$14.1 MILLION by 2027.





NEW Clinical-grade LED Therapy at-home

Non-surgical cosmetic procedure booster recommended by dermatologists The IllumiaSkin PRO uses proven red and blue LED technology for skin rejuvenation to reduce the appearance of fine lines, wrinkles and blemishes.¹⁻⁹



Clinically-proven LED Technology Safe, non-invasive, with multidimensional results

Patented

References: 1. Healthline. LED light therapy for skin: what to know. Available at: healthline.com/health/beauty-skin-care/led-light-therapy (accessed July 2023). 2. American Academy of Dermatology Association. Lasers and lights. How well do they treat acne? Available at: aad.org/public/diseases/acne/derm-treat/lasers-lights (accessed July 2023). **3.** Pei S et al. Indian Dermatol Online J 2015;6(3):145–57. **4.** Goldberg DJ and Russell BA. J Cosmetic Laser Ther 2006;8:71–5. **5.** Lee SY et al. Lasers Surg Med 2007;39(2):180–8. 6. Gold MH et al. J Cosmet Laser Ther 2011;13:308–14. 7. Van Tran V et al. Opt Laser Technol 2021;135:106698. 8. Jagdeo J et al. Lasers Surg Med 2018;50:613–28. 9. Healthline. What is red light therapy, and how does it work? Available at: healthline.com/health/red-light-therapy (accessed July 2023).



Pre-programmed Treatments 10-minute rituals to reset and rejuvenate





WHOLESALE \$269.40

ILLUMIASKIN PRINT ADVERTISING APAN JOURNAL



NEW Clinical-grade LED Therapy at-home



Clinically-proven LED Technology Safe, non-invasive, with multidimensional results

Explore more premium at-home beauty solutions at illumiaskin.com

IllumiaSkin is proudly Australian-owned and run by a dedicated team of #LEDLovers, working to harness the power of science to ignite transformative rituals for people and the planet. Illuminate your skin and support initiatives for positive health, equity, and a sustainable future.

Ready to take your supplementary skin care services to the next level?

Email Debbie Lane at connect@illumiaskin.com or call 0456 696 144.

References: 1. Healthline. LED light therapy for skin: what to know. Available at: healthline.com/health/beauty-skin-care/led-light-therapy (accessed July 2023). 2. American Academy of Dermatology Association. Lasers and lights. How well do they treat acne? Available at: aad.org/public/diseases/acne/derm-treat/lasers-lights (accessed July 2023). 3. Pei S et al. Indian Dermatol Online J 2015;6(3):145-57. 4. Goldberg DJ and Russell BA. J Cosmetic Laser Ther 2006;8:71–5. 5. Lee SY et al. Lasers Surg Med 2007;39(2):180–8. 6. Gold MH et al. J Cosmet Laser Ther 2011;13:308–14. 7. Van Tran V et al. Opt Laser Technol 2021;135:106698. 8. Jagdeo J et al. Lasers Surg Med 2018;50:613–28. 9. Healthline. What is red light therapy, and how does it work? Available at: healthline.com/health/red-light-therapy (accessed July 2023).





Patented Flexible Design Medical-grade silicone, lightweight and portable

The IllumiaSkin PRO uses proven . red and blue LED technology for skin rejuvenation to reduce the appearance of fine lines, wrinkles, redness and blemishes.¹⁻⁹



Pre-programmed Treatments 0-minute rituals to reset and rejuvenate



Product	LED Power	Irradiance	Fluence	Treatment Time	Number of
ILLUMIASKIN PRO	Red (633nm) Near Infrared (830nm) Blue (415nm)	30mW/cm ²	18J cm ²	10 min	270 LEDs (9 red and 90 k bulbs
Omnilux Countour	Red (633nm) Near Infrared (830nm)	30mW/cm ²	18J cm ²	10 min	132 LEDs (66 infrared) end
Omnilux Clear	Red (633nm) Blue (415nm)	30mW/cm ²	18J cm ²	10 min	132 LEDs (66 encased in 6
Omnilux Men	Red (633nm) Near Infrared (830nm & 1072nm)	N/A	N/A	10 min	132 LEDs (66 infrared) end
Dr Dennis Gross Spectralite Faceware Pro	Red (633nm) Near Infrared (880nm) Blue (415nm)	N/A	N/A	3 min	162 LEDs (10
Dr Naomi LED it Glow	Red (630nm) Near Infrared (830nm)	20-30 mw/cm ²	N/A	10 min	156 LEDs (78 red) encased
Aduro	Red (633nm) Near Infrared (830m) Blue (470nm) Green (520nm) Yellow (50% red & 50% green) Orange (84% red & 16% green) Cyan (50% blue & 50% green) Purple (50% blue & 50% red)	30mW/cm ²	18J cm ²	10 min	264 LEDs (6 66 red, 66 b encased in 6
Current Body	Red (633nm) Near Infrared (830nm)	30mW/cm ²	18J cm ²	10 min	132 LEDs (66 infrared) end

f LEDs	Number of Pre-Programmed Treatments	RRP	FDA Cleared	ARTG Registered
90 near-infrared, 90 blue) encased in 90	3	\$449		
56 red and 66 near- ncased in 66 bulbs	1	\$605		
66 red and 66 blue) 66 bulbs	1	\$605		
56 red and 66 near- ncased in 66 bulbs	1	\$605	\checkmark	
00 red and 62 blue)	3	\$709		
78 near-infrared, 78 ed in 78 bulbs	1	\$595		X
66 near-infrared, blue and 66 green) 66 bulbs	11	\$529		
56 red and 66 near- ncased in 66 bulbs	1	\$585		X

RRP	FDA Cleared	ARTG Registered	Photobiological Test
\$449			
\$605			
\$605			
\$605			\checkmark
\$709	\checkmark		N/A
\$595	\checkmark	X	N/A
\$529			
\$585		X	N/A